

HelixCript™ 1st-Strand cDNA Synthesis Kit

Kit Contents

HelixCript™ 1 st -Strand cDNA Synthesis Kit			
Cat. No.	CDNA25 (25rxns)	CDNA50 (50rxns)	CDNA100 (100rxns)
Oligo d(T) ₂₀	0.025ml	0.05ml	0.1ml
Random Hexamers	0.025ml	0.05ml	0.1ml
dNTP Mix (each 10mM)	0.025ml	0.05ml	0.1ml
5x RT Reaction Buffer	0.1ml	0.2ml	0.4ml
0.1M DTT	0.025ml	0.05ml	0.1ml
RNase Inhibitor	0.025ml	0.05ml	0.1ml
<i>Thermo</i> Reverse Transcriptase (200unit/ul)	0.025ml	0.05ml	0.1ml
RNase-free Water	1.0ml	1.0ml	1.5ml
Instructions for Use	1ea	1ea	1ea


Description

HelixCript™ 1st-Strand cDNA Synthesis Kit is a complete system for efficient synthesis of 1st-strand cDNA from total RNA or poly(A)⁺-selected RNA. This kit uses HelixCript™ *Thermo* Reverse Transcriptase, a thermostable and RNaseH negative variant of M-MLV RTase, that is active at a temperature range of 42°C ~ 55°C. At the escalated temperature, the cDNA synthesis is enhanced partially due to the less internal structural formation of the template RNA and increased polymerization activity of RTase. The included RNase inhibitor in this kit protects the template RNA from degradation in the reaction mix and enhances the yield of cDNA. 1st-strand cDNA synthesized with this kit can be directly used as a template in a downstream assay such as PCR or real-time PCR.

Application

Generation of first-strand cDNA for subsequent PCR or Real-time PCR
 cDNA library construction

Storage



Store below -20°C

Shelf life



12 months

Quality Control

Each lot of **HelixCRIPT™ 1st-Strand cDNA Synthesis Kit** was tested against predetermined specifications to ensure consistent product quality.

Protocol

1. Recommended amounts of RNA template and primers for first-strand cDNA synthesis.

- 1) RNAs : Total RNA : 10 ng ~ 5 µg
Poly(A)⁺ RNA : 1 ng ~ 500 ng
- 2) Primers : Oligo-d(T)₂₀ : 0.5 µg or 50 pmoles
Random Hexamer : 50 pmoles
Gene-Specific Primer : 15 ~ 20 pmoles

2. Prepare the "Template/Primer Mixture" by in a clean and sterile microtube.

Components	Volumes (µl)
RNA Template	Xµl
Oligo-d(T) ₂₀ , Random Hexamer, or Gene-Specific Primers	1µl
dNTP Mix (each 10 mM)	1µl
RNase-free Water	Yµl
Total	13µl

3. Incubate at 65°C for 5 minutes and immediately place on ice.

4. Prepare the reaction mixture by adding the following components in the indicated order.

Adding order	Components	Volume (µl)
1	Template/Primer Mixture (from step 3)	13µl
2	5x RT Reaction Buffer	4µl
3	<i>Thermo</i> Reverse Transcriptase (200unit/ul)	1µl
4	RNase inhibitor	1µl
5	0.1M DTT	1µl
Total		20µl

* (Caution) If you do not follow the indicated order, precipitations or loss-of-activity of enzymes will be occurred.

5. Incubate the reaction mixture at 50°C for 30 ~ 50 minutes for gene-specific primer.

When oligo-d(T) or random hexamer is used in reaction, perform 10 min at 42°C, followed by 30 ~ 50 min at 50°C.

6. Inactivate the reaction by heating at 70 °C for 10 min.
7. Synthesized cDNA is immediately used for PCR or store at -20 °C .

Products

Cat. No.	Products	Size
CDNA25	HelixCript™ 1 st -Strand cDNA Synthesis Kit [All-in-One]	25rxns
CDNA50	HelixCript™ 1 st -Strand cDNA Synthesis Kit [All-in-One]	50rxns
CDNA100	HelixCript™ 1 st -Strand cDNA Synthesis Kit [All-in-One]	100rxns