

# HelixAmp™ Hot-Start Fidelity Polymerase

## Kit Contents

HelixAmp™ Hot-Start Fidelity Polymerase	
Cat. No.	HSF250 (250units)
Hot-Start Fidelity Polymerase (2.5unit/μl)	0.1ml
5x HSF Buffer	1.0ml x 2ea
5x TuneUp™ Solution	1.0ml
Instructions for Use	1ea

## Description

**HelixAmp™ Hot-Start Fidelity Polymerase** is a hot-start formulation of a modified *Pfu* DNA Polymerase, meticulously optimized for robust, high-fidelity, and specific amplification of DNA fragments up to 5kb in length. This high-fidelity enzyme demonstrates a 2.5-fold increase in fidelity compared to conventional *Pfu* DNA polymerase and 100-fold higher than *Taq* DNA polymerase in accuracy. The provided 5x HSF buffer (a reaction buffer) includes all essential PCR components, such as buffering reagents, magnesium, salts, and dNTPs. Additionally, the supplied TuneUp™ Solution is for the challenging amplifications of templates with high G+C content or structural complexities, such as repeat sequences.

## Application

Gene cloning with blunt ends  
 High Fidelity/Specific PCR  
 NGS template generation

## Storage buffer

50mM Tris-HCl (pH 8.0), 100mM KCl,  
 0.1mM EDTA, 1mM DTT, stabilizers,  
 50% Glycerol

### Storage



Store below -20°C

### Shelf life



12 months

### Concentration

2.5 unit/μl

## Quality Control

Each lot of **HelixAmp™ Hot-Start Fidelity Polymerase** was tested against predetermined specifications to ensure consistent product quality.

## Protocol

### 1. Recommended amount of template DNA

Human genomic DNA : 10 ~ 500ng / *E.coli* genomic DNA : 100 ~ 200ng

Purified plasmid DNA : 10pg ~ 10ng / cDNA : 25 ~ 750ng

### 2. Mix following components in a PCR tube

Components	Volumes (μl)
Template	X μl
5x HSF Buffer	10 μl
Forward Primer (10uM)	1 ~ 2 μl
Reverse Primer (10uM)	1 ~ 2 μl
5x TuneUp™ Solution [optional] <sup>1)</sup>	0 ~ 20 μl
Hot-Start Fidelity Polymerase (2.5unit/ul)	0.5 μl
RNase-free Water	to 50 μl

<sup>1)</sup> The application of TuneUp™ Solution is beneficial for amplifying challenging templates, including those with high G+C content and repeat sequence regions. To integrate TuneUp™ Solution into the PCR reaction mixture, it is added at a final concentration ranging from 0.5x to 2x.

### 3. PCR condition

Temperature & time	Cycles
95°C, 5 min	x 1
95°C, 20 sec	} x 25 ~ 40
Annealing Temp., 30 sec	
72°C, 1min/kb (Expected size of product)	
72°C, 3 min	x 1

## Products

Cat. No.	Products	Size
HSF250	HelixAmp™ Hot-Start Fidelity Polymerase	250units