

HelixAmp™ Fast RT-LAMP Lyo-Cake (V3)

Kit Contents

HelixAmp™ Fast RT-LAMP Lyo-Cake (V3)		
Cat. No.	LFRLP3-C96 (96rxns)	LFRLP3-C480 (480rxns)
Rehydration Buffer	1.0ml x 2ea	1.0ml x 10ea
Fast RT-LAMP Lyo-Cake (V3)	8-strip x 12ea	8-strip x 60ea
8-Cap Strip	12ea	60ea
Instructions for Use	1ea	1ea

Description

HelixAmp™ Fast RT-LAMP Lyo-Cake (V3), a lyophilized RT-LAMP reagent, provides a one-step solution for loop-mediated isothermal amplification (LAMP) of RNA targets. Especially this kit suppress non-specific product formation in isothermal amplification. Through the lyophilization method, moisture is removed from the product to facilitate storing or transporting products at room temperature. The Lyophilized RT-LAMP premix can be rehydrated simply by adding the Rehydration Buffer along with primers and templates.

HelixAmp™ Fast RT-LAMP Lyo-Cake (V3), contains engineered *Bst* DNA Polymerase, thermo-stable reverse transcriptase (RTase), RNase inhibitor, dNTPs, Mg²⁺, salts, and stabilizing agents. This kit allows the fast amplification of the target within 30 minutes with an appropriate LAMP primer set. The thermo-stable RTase is fully active at a relatively high temperature (60°C), making the one-step RT-LAMP possible at a constant reaction temperature.

Application

Loop-mediated Isothermal Amplification (LAMP) of RNA

Storage

- Store at room temperature for a maximum of 6 months.
- For extended storage, maintain the product at or below 4°C for up to 12 months.
- Please protect it from direct light exposure.
- After opening the pouch, we suggest promptly resealing it and storing it in a freezer until the next use.

Quality Control

Each lot of **HelixAmp™ Fast RT-LAMP Lyo-Cake (V3)** was tested against predetermined specifications to ensure consistent product quality.

Protocol

1. Reaction Mixture

- (1) Remove the seal of the lyophilized cake tube.
- (2) Add the reaction components to the lyophilized RT-LAMP cake tube as in the following table for the selected analysis method (end-point or real-time assay) :

Components	For end-point	For real-time (Intercalating dye)
RNA Template	1 ~ 5 μ l	1 ~ 5 μ l
10x LAMP Primer Mix ¹⁾	2 μ l	2 μ l
Intercalating dye ²⁾	-	x μ l
Rehydration Buffer	Adjust to final 20 μ l. ※ Use at least 10 μ l of Rehydration Buffer.	

- ¹⁾ The following table lists the recommended concentrations for each primer in the 10x LAMP Primer Mix. If you encounter low efficiency or non-specific amplification, adjust the primer concentrations or design a new set of primers for the target sequence.

10x LAMP Primer Mix	
LAMP primers	Primer concentration.
FIP	16 μ M
BIP	16 μ M
F3	2 μ M
B3	2 μ M
LF	8 μ M
LB	8 μ M

- ²⁾ It is recommended to use SYBR Green I (final concentration of 0.1-0.2x), EvaGreen dye (0.1-0.3x), or SYTO 9 (0.1-0.5x). Please note that intercalating fluorescent dyes are not included in this kit.
- (3) Seal the tube with a cap strip. Gently tap the reaction mix and centrifuge briefly.
- (4) Perform the LAMP reaction according to the below Reaction Condition.

2. Reaction Condition

For end-point assay: Incubate at 60°C for 30 minutes. Time can be extended as necessary for low-copy targets or challenging sample types.

For real-time assay (fluorescent dye): Use a real-time PCR machine or an isothermal amplification instrument to run the assay. Set the instrument to a constant incubation temperature of 60°C. Measure the fluorescence intensity every 1 min for 30 minutes. The reaction time can be increased as necessary for low copy targets or challenging sample types.

Products

	Products	Cat. No.	Size
Cake type	HelixAmp™ FastLAMP Lyo-Cake (V3)	LFLP3-C96	96rxns
		LFLP3-C480	480rxns
	HelixAmp™ Fast RT-LAMP Lyo-Cake (V3)	LFRLP3-C96	96rxns
		LFRLP3-C480	480rxns
Bead type	HelixAmp™ FastLAMP Lyo-Dot (V3)	LFLP3-B96	96rxns
	HelixAmp™ Fast RT-LAMP Lyo-Dot (V3)	LFRLP3-B96	96rxns