

HelixCRIPT™ One-Step RT-PCR Kit [*Hot-Taq*]

Kit Contents

HelixCRIPT™ One-Step RT-PCR Kit [<i>Hot-Taq</i>]		
Cat. No.	ORTHT100 (100 rxns)	ORTHT500 (500 rxns)
Enzyme Mix [<i>Hot-Taq</i>]	0.2 ml	0.2 ml x 5ea
2x Reaction Mix [<i>Hot-Taq</i>] (containing dNTP mix, MgCl ₂)	1.25 ml x 2ea	1.25 ml x 10ea
Instruction for Use	1ea	1ea

* Store at -20°C

Description

HelixCRIPT™ One-Step RT-PCR Kit [*Hot-Taq*] is designed for sensitive amplification of the target gene in one-tube reaction from total transcripts. A reverse transcriptase, HelixCRIPT™ *Thermo* Reverse Transcriptase, and a HelixAmp™ *Hot-Taq* polymerase are supplied as an enzyme mixture.

HelixCRIPT™ One-Step RT-PCR Kit [*Hot-Taq*] : High specific and sensitive RT-PCR

One-step RT-PCR system provide the several advantages.

- Synthesis of cDNA and PCR amplification corresponding to target gene in one-tube reaction
- Obtain the reproductive data in the repetitive experiment(s)
- Can save the time and cost for preparation of RT-PCR
- Amplification of low-copy transcripts by RT-PCR

Application

Detection of target gene transcript from RNA

Semi-quantitative, quantitative analysis of RNA transcription level

Quality control assay data

Functional analysis

The activity for cDNA synthesis and PCR amplification of target gene transcript using HelixCRIPT™ One-Step RT-PCR Kit was evaluated by Limit-of Detection (LOD) assay and long range of gene constitutively-expressed in human total transcripts.

Quality authorized by Yountaek Go



NanoHelix Co., Ltd.

43-15, Techno 5-ro, Yuseong-Gu, Daejeon, 34014, South Korea. TEL : 82-42-867-9055, FAX : 82-42-867-9057

E-mail : info@nanohelix.net < www.nanohelix.net www.nanohelix.co.kr/KOR >

Protocol

1. Program the thermal cycler as follows in order to synthesize cDNA using HelixCRIPT™ One-Step RT-PCR Kit with HelixAmp™ *Hot-Taq*.

Step	Condition		Cycle(s)
cDNA Synthesis	42 ~ 55°C for 30 ~ 50 min		1
Pre-denaturation	95°C for 12 ~ 15 min		1
PCR Amplification	Denaturation	95°C for 20 sec	30 ~ 40
	Annealing	¹⁾ AT°C for 40 sec	
	Extension	72°C for 1 min/kb Collect the fluorescence data	
Post extension	72°C for 5 min		1

¹⁾ AT, annealing temperature of primers used

$$\text{Annealing Temperature} = T_m - (6 \sim 8^\circ\text{C})$$

$$\text{Where, } T_m (\text{Melting Temp.}) = [4^\circ\text{C} \times (\text{G} + \text{C})] + [2^\circ\text{C} \times (\text{A} + \text{T})]$$

2. Add the following components into 0.2 or 0.5 ml micro-tube.

Components	Volumes
2x Reaction Mix [<i>Hot-Taq</i>]	25 μl
RNA Template (1 ng ~ 5 μg)	X μl
Forward primer (10 pmoles/ μl)	2 μl
Reverse primer (10 pmoles/ μl)	2 μl
Enzyme Mix [<i>Hot-Taq</i>]	2 μl
RNase-free water	to 50 μl

※ RNAs : Total RNA : 10 ng ~ 5 μg , Poly(A)⁺ RNA : 1 ng ~ 500 ng

3. Gently mix and immediately centrifuge the reaction mix.
4. Perform the one-step RT-PCR.

Products

Cat. No.	Products	Size
ORTHT100	HelixCRIPT™ One-Step RT-PCR Kit [<i>Hot-Taq</i>]	100 rxns
ORTHT500	HelixCRIPT™ One-Step RT-PCR Kit [<i>Hot-Taq</i>]	500 rxns

NanoHelix Co., Ltd.

43-15, Techno 5-ro, Yuseong-Gu, Daejeon, 34014, South Korea. TEL : 82-42-867-9055, FAX : 82-42-867-9057

E-mail : info@nanohelix.net < www.nanohelix.net www.nanohelix.co.kr/KOR >