

HelixAmp™ Premium-Pfu Polymerase

Kit Contents

| HelixAmp™ Premium-Pfu Polymerase | | | |
|----------------------------------|---------------------|---------------------|------------------------|
| Cat. No. | PP250 (250units) | PP500 (500units) | PP2500 (2,500units) |
| Premium-Pfu (2.5unit/μl) | 0.1ml | 0.2ml | 0.2ml x 5ea |
| 5x Reaction Mix [PP] | 1ml x 2ea | 1ml x 4ea | 1ml x 20ea |
| N-Solution™ | 1ml | 1ml x 2ea | 1ml x 10ea |
| 6x Loading Dye | 0.5ml | 1ml | 1ml x 5ea |
| Blue Box | - | - | 1ea |
| Instruction for Use | 1ea | 1ea | 1ea |

※ The 5x Reaction Mix [PP] contains buffers, dNTPs, Mg²⁺, and BSA in appropriate concentrations.

Description

HelixAmp™ Premium-Pfu Polymerase is an editorial enzyme improved from HelixAmp™ Speed-Pfu Polymerase. HelixAmp™ Premium-Pfu Polymerase shows faster (4 ~ 6 times) polymerization reaction and higher fidelity (3 ~ 5 times) than any other high-fidelity enzymes. NanoHelix's "PMT (polymerase modulator on temperature) technology" is applied in the buffer system, which is effective to reduce primer-dimer formation and non-specific amplification during the PCR. HelixAmp™ Premium-Pfu Polymerase is most suitable to faithful amplification of relatively long-ranged target for cloning etc. Due to its high speed, fast PCR with this enzyme could be completed in 30 min for the reliable amplification of less than 1 kb size target DNA. For the maximum performance of PCR reactions high-quality dNTP mixture is supplied. In case of PCR amplification of target DNA with high G+C content or structural problem, such as repeat sequence, the application of N-Solution™ improves the specificity and productivity of the reactions.


Application

High fidelity PCR
 Long range PCR
 Site-directed mutagenesis
 Blunt end PCR cloning

Storage buffer

50mM Tris-HCl (pH 8.0), 0.1mM EDTA,
 1mM DTT, stabilizers, 50% Glycerol

Storage



Shelf life



Concentration

2.5unit/μl

NanoHelix Co., Ltd.

A-dong and B-dong, 43-15, Techno 5-ro, Yuseong-Gu, Daejeon, 34014, South Korea. TEL : 82-42-867-9055, FAX : 82-42-867-9057

E-mail : info@nanohelix.net <www.nanohelix.net www.nanohelix.co.kr/KOR>

Quality Control Assay

Contamination Assay

HelixAmp™ Premium-Pfu Polymerase was passed from quality control assay for contamination of bacterial host DNA.

Functional assay

HelixAmp™ Premium-Pfu Polymerase was functionally tested for PCR amplifications using the various primer sets (0.5 kb ~ 17kb) from human *beta-globin* gene.

Quality authorized by Yountaek Go



Protocol

1. Recommended amount of template DNA.

Human genomic DNA : 10 ~ 100ng

Bacterial genomic DNA : 5 ~ 50ng

Purified plasmid or phage DNA : 1 ~ 5ng

2. Mix following components in a PCR tube.

| Components | Reaction Volumes (μl) | |
|------------------------------|-----------------------|-----------|
| | 25 | 50 |
| Template | X μl | X μl |
| Forward Primer (10pmoles/μl) | 1μl | 2μl |
| Reverse Primer (10pmoles/μl) | 1μl | 2μl |
| 5x Reaction Mix [PP] | 5μl | 10μl |
| N-Solution™ [optional]※ | 0~2.5μl | 0~5μl |
| Premium-Pfu | 1.25units | 1.25units |
| RNase-free Water | to 25μl | to 50μl |

※ **N-Solution™** is an additive altering the binding behavior of primer and template and can help the amplification that do not work well under standard PCR condition. Especially, **N-Solution™** can be used for the amplification of problematic template, such as high G+C content and repeat sequence regions. The optimal concentrations of **N-Solution™** are vary upon the primer-template sets and should be set by adding into the PCR reaction mixture from 2 to 10% volume. Most of the PCR reactions are not required the **N-Solution™** and we recommend to use the **N-Solution™** only in case of the PCR amplification is not works well or too much non-specific products are observed.

NanoHelix Co., Ltd.

A-dong and B-dong, 43-15, Techno 5-ro, Yuseong-Gu, Daejeon, 34014, South Korea. TEL : 82-42-867-9055, FAX : 82-42-867-9057

E-mail : info@nanohelix.net <www.nanohelix.net www.nanohelix.co.kr/KOR>

3. PCR condition.

| Temperature & time | Cycles |
|--|-------------|
| 95°C, 2 min | x 1 |
| 95°C, 20 sec | } x 25 ~ 40 |
| Annealing Temp. ¹⁾ , 40 sec | |
| 72°C, 30 sec/kb (Expected size of product) | |
| 72°C, 5 min | x 1 |

¹⁾ Annealing Temp. = $T_m - (4 \sim 6^\circ\text{C})$

T_m (Melting Temp.) = $[4^\circ\text{C} \times (\text{number of G} + \text{C})] + [2^\circ\text{C} \times (\text{number of A} + \text{T})]$

※ If there are non-specific PCR amplifications, use an annealing temperature of 2~5°C higher than the calculated melting temperature.

Products

| Cat. No. | Products | Size |
|----------|----------------------------------|------------|
| PP250 | HelixAmp™ Premium-Pfu Polymerase | 250units |
| PP500 | HelixAmp™ Premium-Pfu Polymerase | 500units |
| PP2500 | HelixAmp™ Premium-Pfu Polymerase | 2,500units |

NanoHelix Co., Ltd.

A-dong and B-dong, 43-15, Techno 5-ro, Yuseong-Gu, Daejeon, 34014, South Korea. TEL : 82-42-867-9055, FAX : 82-42-867-9057

E-mail : info@nanohelix.net <www.nanohelix.net www.nanohelix.co.kr/KOR>