

## RealHelix™ qPCR Kit [Green]

### Kit Contents

RealHelix™ qPCR Kit [Green]		
Cat. No.	QP2-S200 (200rxns)	QP2-S500 (500rxns)
2x Premix [Green]	1ml x 2ea	1ml x 5ea
ROX Dye (25µM)	0.2ml	0.5ml
Instructions for Use	1ea	1ea


### Description

**RealHelix™ qPCR Kit [Green]** is designed to perform a rapid real-time quantification of target DNA using a double strand DNA-binding green fluorescent dye. The convenient 2x premix contains hot-start PCR enzyme, green fluorescent dye, dNTPs, buffers, Mg<sup>2+</sup>, and stabilizing agent. The hot-start PCR enzyme provides high specific amplification of target DNA and minimizes the side products such as primer dimers. The separately supplied ROX Dye could be used as a passive reference dye in real-time PCR instruments that are compatible with the evaluation of the ROX signal.

### Application

Quantification of target DNA

#### Storage



In the dark at -20°C.

#### Shelf life



12 months

## Quality Control

By Nanohelix's ISO 13485-certified quality management system, each lot of **RealHelix™ qPCR Kit [Green]** was tested against predetermined specifications to ensure consistent product quality.

## Protocol

1. Program a real-time PCR instrument according to the recommendations below. The excitation and emission maxima of the included green fluorescence dye are at 494 nm and 521 nm, respectively.

### <2-step cycling method>

Step	Condition		Cycle(s)
Enzyme Activation		95°C for 15 min	1
PCR Amplification	Denaturation	95°C for 20 sec	40
	Annealing & Extension	60°C for 40 sec <b>Collect the fluorescence data</b>	

### <3-step cycling method>

Step	Condition		Cycle(s)
Enzyme Activation		95°C for 15 min	1
PCR Amplification	Denaturation	95°C for 20 sec	40
	Annealing	<sup>1)</sup> AT°C for 30 sec	
	Extension	72°C for 1 min/kb <b>Collect the fluorescence data</b>	

<sup>1)</sup> **AT**, annealing temperature of primers used

Annealing Temperature =  $T_m - (6 \sim 8^\circ\text{C})$

Where,  $T_m$  (Melting Temp.) =  $[4^\circ\text{C} \times (\text{G} + \text{C})] + [2^\circ\text{C} \times (\text{A} + \text{T})]$

☞ For melting curve analysis of products, refer the instrument documentation

## 2. Add following components for a 20 $\mu$ l volume reaction.

Components	Volumes
DNA Template	X $\mu$ l
2x Premix [Green]	10 $\mu$ l
Forward primers (10 $\mu$ M)	0.5 $\mu$ l (final 0.25 $\mu$ M) <sup>1)</sup>
Reverse primers (10 $\mu$ M)	0.5 $\mu$ l (final 0.25 $\mu$ M) <sup>1)</sup>
ROX Dye (25 $\mu$ M)	Optional <sup>2)</sup>
RNase-free Water	Adjust to final 20 $\mu$ l

<sup>1)</sup> The amount of each primer should be adjusted according to the efficiencies of target amplification.

<sup>2)</sup> Use the recommended amount or concentration of ROX Passive Reference Dye depending on the instrument.

## 3. Gently mix and immediately centrifuge the reaction mix.

## 4. Perform the Real-time PCR.

## Products

Cat. No.	Products	Size
QP2-S200	RealHelix™ qPCR Kit [Green]	200rxns
QP2-S500	RealHelix™ qPCR Kit [Green]	500rxns