

RealHelix™ Superplex qPCR Kit [Probe]

Kit Contents

RealHelix™ Superplex qPCR Kit [Probe]		
Cat. No.	SUQ-P200 (200rxns)	SUQ-P500 (500rxns)
Superplex 2x Premix [Probe]	1ml x 2ea	1ml x 5ea
ROX Dye (25µM)	0.2ml	0.5ml
Instructions for Use	1ea	1ea


Description

RealHelix™ Superplex qPCR Kit [Probe] is designed for multiple target qPCR detections using dual-labeled probes in real time PCR instruments. The kit consists of convenient 2x concentrated premix solution and ROX Dye. The user convenient 2x premix contains all required components including a hot-start Taq polymerase, heat-labile Uracil-DNA glycosylase, dUTP, dNTPs, buffer, magnesium ion and stabilizing agents. The hot-start Taq polymerase provides high specific amplification of target DNA and minimizes the side products such as primer dimers. This UDG/dUTP system prevents the carryover contamination of PCR products from previous reactions. Included ROX Dye could be used as a passive reference dye.

Application

Quantification of target DNA sample by real time PCR

Storage



Store below -20°C

※ ROX Dye should be stored in the dark.

Shelf life



12 months

Quality Control

By Nanohelix's ISO 13485-certified quality management system, each lot of **RealHelix™ Superplex qPCR Kit [Probe]** was tested against predetermined specifications to ensure consistent product quality.

Protocol

1. Program a real-time PCR instrument according to the recommendations below.

Step	Condition	Cycle(s)	
[Optional] UDG Reaction*	20 ~ 25°C for 5 min	1	
Enzyme activation	95°C for 15 min	1	
PCR Amplification	Denaturation	40	
	Annealing**		AT°C for 20~30 sec
	Extension		72°C for 1 min/kb
Collect the fluorescence data			

* The UDG reaction step is not essential. The UDG will efficiently remove carryover contaminant DNA during sample setup and cyclor ramping.

** Annealing Temperature = $T_m - (6 \sim 8^\circ\text{C})$

Where, T_m (Melting Temp.) = $[4^\circ\text{C} \times (\text{G} + \text{C})] + [2^\circ\text{C} \times (\text{A} + \text{T})]$

2. Add following components for a single 20µl reaction volume. *

Components	Volumes
DNA Template (1ng ~ 1µg)	x µl
Superplex 2x Premix [Probe]	10µl
Forward primers (10µM)	0.2 ~ 0.5µl ¹⁾
Reverse primers (10µM)	0.2 ~ 0.5µl ¹⁾
Probes (10µM)	0.1 ~ 0.5µl
ROX Dye (25µM)	Optional ²⁾
RNase-free Water	Adjust to final 20µl

* For multiple reactions, prepare a master mix by adding the required volumes of each above components (except the DNA template) and dispense appropriate volumes into PCR tubes or plates.

- 1) The amount of each primer should be adjusted according to the efficiencies of target amplifications.
- 2) Use the recommended amount or concentration of ROX Dye (Passive Reference) depending on the instrument.

3. Gently mix and briefly centrifuge the reaction mix.

4. Perform the Real-time PCR.

Products

Cat. No.	Products	Size
SUQ-P200	RealHelix™ Superplex qPCR Kit [Probe]	200rxns
SUQ-P500	RealHelix™ Superplex qPCR Kit [Probe]	500rxns